


New nucleic acids extraction technologies from FFPE fixed tissues




» Get rid of the deparaffinization step

- No deparaffinization step, no overnight lysis
- From your fixed samples to your qPCR results within the day



» Yield and quality

- High yields, from small to long fragments of RNA
- High quality RNA from FFPE samples, slides or microdissected
- Up to 5 slides (10 µm thick, 1 to 2 cm²) per samples, 2 to 6 µg of RNA per slide




» Flexibility with the AmpTec technology

- On-column or magnetic extractions
- Total RNA extraction
- Separate miRNA from mRNA

» Amplify you RNA if needed

Unique amplification technology for your FFPE samples



ExcilONE

Service for laser microdissection and nanogenomic



6, Rue Blaise Pascal - Parc Euclide - Bâtiment A - 78990 Elancourt - France
Tél. : +33 (0)1 34 86 77 01 • Fax : +33 (0)1 34 86 77 53
Email : service@excilone.com • Internet : www.excilone.com